Sys 2128227

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









RUA TUPI N° 535 S. PAULO

TELS 67-7806-05-0912 • 676-6722

# Laser Doppler Anemometry Equipment Catalog

## Contents Page

Using LDA in Flow Measurements	2
55X Modular LDA Optics	3
55L10 Photomultiplier	6
55L11 Diode Detector	7
55L67 Optical Unit	7
55L17 Seeding Generator	7
LDA Signal Processing	7
55N20 Doppler Frequency Tracker	8
55L90a LDA Counter Processor	10
55N10 Frequency Shifter	12
57G20 Buffer Interface	13
Educational LDA System	14
Low Velocity LDA System	14
Coaxial Laser System	14
Backscatter LDA System	15
Two-Color LDA System	15
LDA Computer System	16

*Handwritten notes:*  
14 580  
3/5

# Using LDA in Flow Measurements

## Introduction

In many cases, a flow can only be measured correctly if a non-contact sensor is used. For example, a flow pattern may be of so small dimensions that even the smallest of probes would affect the parameters under investigation, and a flow of a chemically active medium would damage a physical probe. Some flows even make the physical presence of a probe impossible; cases in point are flows around propellers, flows in turbines, etc.

In such flows Laser Doppler Anemometry (LDA) may be the only possible means of measuring local flow velocity, and in many other flows LDA may be the most practical method of investigating the dynamic properties of flows.


Fig. 1. Vectors in LDA

## LDA Theory

The Laser Doppler Anemometer uses the Doppler shift (Christian Doppler, 1842) of light scattered by moving particles to determine particle velocity and thus find the fluid flow velocity. The general equation expressing the Doppler shift ( $f_D$ ) in the frequency of the scattered light as a function of particle velocity  $\vec{V}$  is:

$$f_D = f_s - f_i = \frac{1}{\lambda} \cdot \vec{V} \cdot (\hat{e}_s - \hat{e}_i)$$

$f_s$ , frequency of scattered light  
 $f_i$ , frequency of incident light  
 $\hat{e}_s, \hat{e}_i$ , unit vectors of scattered and incident respectively and  
 $\lambda$ , wavelength of incident light

The use of a monochromatic coherent light source such as a laser beam makes it possible to determine the values of  $f_i$  and  $\lambda$ .  $\hat{e}_s, \hat{e}_i$  depend on the geometry of the system. The velocity component  $V_x$  measured by LDA will be in the plane formed by  $\hat{e}_i$  and  $\hat{e}_s$  and perpendicular to the line bisecting the angle formed by the two unit vectors (Fig. 1).

## LDA Modes

Three modes of operation - the Reference Beam Mode, the Differential Doppler Mode, and the Dual Scatter Mode - have been used in Laser Doppler Anemometry, but only the first two have found general acceptance.


Fig. 2. Reference Beam Mode

**Reference Beam Mode.** This mode was used by YEH and CUMMINGS in their original article in Applied Physics, 1964, in which LDA was first presented. In the Reference Beam Mode light scattered from one laser beam mixes in the photodetector with light from a reference beam. In practice this is done by a coaxial superposition by means of beam splitters and mirrors of the scattered beam and the reference beam at the photodetector (Fig. 2). The reference beam impinges constantly on the photodetector, and best S/N ratio will be obtained at high particle densities. Typical signal characteristics are: Continuous modulation with coherent parts separated only by amplitude dropouts at phase shifts, good to moderate S/N ratio at high particle densities (water flows), and poor S/N ratio at low particle densities (air flows). Since the equation for  $f_D$  includes the unit vector  $\hat{e}_s$  for scattered light, the photodetector must have a small collecting aperture and be placed coaxially with the reference beam.

## Differential Doppler Mode

(Fig. 3). This mode uses two beams of equal intensity as scattering beams. For any direction of detection the light picked up by the photodetector will consist of the superposition of light scattered from either beam.  $f_D$  can be found as the difference between the two scattering frequencies  $f_{s1}$  and  $f_{s2}$ , and it can be seen that  $f_D$  is independent of direction of scattering  $\hat{e}_s$ . Because of the independence of the positioning of the receiving optics, the Differential Doppler Mode can be used both in forward scatter and in backscatter, making it possible to adapt the LDA system to different configurations. It can also be used with a large receiving aperture, enabling measurement on small single particles.


Fig. 3. Differential Doppler Mode

Typical signal characteristics are: Good S/N ratio from single particles with low noise level between bursts (noise proportional to light level on photodetector), good-to-moderate S/N ratio on flows of high particle density (air or water), poor S/N ratio on flows having high density of very small particles (typically, polluted water flow).

**Fringe Model.** The Differential Doppler Mode can be explained more directly by means of a


Fig. 4. Fringe Model

**Fringe Model (Fig. 4).** At the beam crossing, the two laser beams can be described as two beams of monochromatic coherent light, each with plane and parallel wave fronts. In the intersection, alternately constructive and destructive superposition of the two beams will form a fringe pattern of planes having a strong electromagnetic field separated by planes having a weak electromagnetic field. Distance between fringes,  $d_f$ , can be calculated from the angle of intersection  $\theta$  and wavelength  $\lambda_c$  of the laser light.

$$d_f = \frac{\lambda_c}{2 \sin \theta/2}$$


A particle moving across the fringe pattern with a velocity component  $V_x$ , perpendicular to the fringe planes will scatter light at a frequency  $f_D$ .

$$f_D = \frac{V_x}{d_f} = V_x \cdot \frac{\theta}{2 \sin \theta/2}$$

# 55X Modular LDA Optics

## Features:

- ★ Fully modular optical system
- ★ Coaxial design
- ★ Measured velocity vector plane can be rotated through 360°
- ★ Hard anti-reflex coating on all optical surfaces
- ★ Optical filters and coating rated for high power laser
- ★ Optics corrected for spherical and chromatic aberrations
- ★ Ball-bearing supports featuring graticules
- ★ Beam expansion and full-aperture receiving optics for higher signal-to-noise ratio
- ★ Adjustable receiving optics for precise focusing
- ★ Only one frequency shift module and one backscatter module for one- and two-dimensional systems
- ★ Mounting facilities for both HeNe-lasers and Argon-Ion lasers

**Introduction**  
 The DISA 55X Modular LDA Optics constitutes a third generation in the development and design of laser Doppler optics. The overall optical layout is optimized for the utmost performance based on intensive research and developing experience from all types of LDA applications. Great care has been taken to develop the individual optical components (beam splitters, color separators, filters, etc.) utilizing advanced multilayer technique.


Some typical LDA Systems for setting up the 55X Optics for your application are shown on pages 14 and 15, with information on mode of operation. To start it, select the optical bench 55X19, 55X41 or 55X42 attaching the laser head to be

used. (The benches comprise mounting parts). Select one or two supports (55X23) for mounting the various modules, such as beam splitter, Bragg Cell, or beam translator. The modules feature standard lengths, inside screws, and guide pins for perfect mounting. Complete your LDA optics with a cover at

the laser end, and a mounting ring with the actual front lens facing your investigation setup. In forward scatter modes the photomultiplier can be placed on the tripod (or clamped to your investigation setup), and in the backscatter mode be mounted together with the backscatter module.

*The 55X Modular LDA Optics comprises a complete line of optical components for all modes of operation occurring in practice.*

*Two-Color Backscatter System*

# X Modular LDA Optics


**55X08 PM Section**  
This photomultiplier section uses the PM tube and has a padded mounting ring for the 55X36/37/38 Interference Filters. In a two-color setup two sections with Interference Filters are mounted to the 55X35 or Separator, which in turn mounts to the 55X34 PM Optics. Photomultiplier Tube RCA 4526  
x. Cathode Voltage 2000 V, negative polarity  
x. Anode Current 0.1 mA  
Spectral Response  
Ext. S20, quantum efficiency 14% (red) 20% (green/blue)  
Sensitivity 4400 A/W at 550 nm  
Bandwidth 0 to 120 MHz, 50 Ω  
Dimensions 120 mm length, 52 mm dia  
Weight 0.5 kg

**55X19 Laser Adapter**  
The Laser Adapter comprises a holder for a coaxial laser head and a mounting bench with adjustable feet for the complete LDA Optics.

Coaxial Laser Diameters  
50 mm or 44.5 mm  
Dimensions  
Adaptor tube: 165 mm length, 112/60 mm dia  
Mounting Bench: 360 mm length, 108 mm width  
Weight  
Adaptor tube: 0.5 kg  
Mounting bench: 1.8 kg

**55X22 Beam Waist Adjuster**  
The Beam Waist Adjuster is used to optimize the fringe pattern quality in the measuring volume. It is important to adjust the laser beam divergence in a system having a front lens with long focal length. The 55X22 Adjuster mounts in the 55X20/21 Cover. Adjustable Focal Length -0.8 m to +0.5 m, beam expansion x 1.1  
Dimensions 50 mm length, 43 mm dia

**55X24 Beam Splitter**  
This beam splitter module has a polarizing beam splitter prism, a polarisator rotator and a prism. The module has variable beam splitter ratio and is used in reference beam mode as well as in differential Doppler mode with photomultiplier. When operating in reference beam mode together with the 55L11 Diode Detector the polarisator rotator should be out of function (pull the small rod.)  
Beam Splitter  
polarization prism  
633 nm  
Polarization Rotator  
90°, 633 nm  
Beam Separation  
One beam coaxial, one displaced 30 mm  
Dimensions  
35 mm width, 112 mm dia  
Weight 0.7 kg

**55X12 Beam Expander**  
The Beam Expander can be used in all 55X systems with variable beam separation (55X32 Beam Separator). The Beam Expander is essential for applications in which the size of the measuring volume is critical. By expanding the laser beam the diameter of the measuring volume can be halved, resulting in quadrupling of the light intensity.

Beam Expander Ratio  
95:50  
Beam Expander Optics  
Achromatic with hard AR-coating. R<0.3% per surface  
Dimensions  
170 mm length, 112 mm dia  
Weight 1.0 kg


**55X20/21 Cover and Retarder**  
This set of Cover and Retarder is available for Argon-Ion laser 55X20 and for HeNe-laser 55X21. The λ/4 retarder is mounted on the laser head, and the cover with retarder polarized laser beam will after passing these components have its polarization plane following the 55X Optics when this is rotated.

**55X20 Cover and Retarder**  
Wavelength of light: 488 to 515 nm

**55X21 Cover and Retarder**  
Polarization can be rotated relative to optics  
Wavelength of light: 633 nm  
Thread of Retarder λ/4 type  
UNF 1" x 32 (laser head standard)


Dimensions  
Cover: 112 mm dia  
Retarder: 37 mm dia  
Weight 0.3 kg

**55X23 Support**  
One Support is used for mounting small 55X systems and two supports secure perfect stability of large optical systems. The optical benches have mounting holes for the supports.  
Rotatability  
360°, graticule ± 90°, ball bearing  
Optical Axis  
110 mm over bench  
Dimensions  
170 mm height, 126 mm, 17.5 mm thick  
Weight 1.1 kg

**55X25 Beam Splitter**  
The 55X25 Beam Splitter has a neutral beam splitter and a prism, and is mainly used in two color setups.  
Beam Splitter  
Reflecting prism type  
Beam separation  
One coaxial beam, one beam separated 30 mm  
Dimensions  
35 mm width, 112 mm dia  
Weight 0.7 kg

#### 55X26 Beam Splitter

The 55X26 Beam Splitter is of a neutral type and splits an incoming coaxial beam into two beams parallel to the axis. This module is mainly used in 2-channel reference beam mode, and for one component measurement without frequency shift.

#### Beam Splitter

#### Neutral prism type

#### Beam Separation

Each beam separated 30 mm from optical axis, in the same plane

#### Dimensions

35 mm width, 112 mm dia

#### Weight

0.7 kg

#### 55X28 Beam Displacer

The 55X28 Beam Displacer is used together with the 55X24 Beam Splitter in one-component setups (also with frequency shift) Beam Displacement

Incoming coaxial beam is displaced 30 mm by prisms

#### Dimensions

35 mm thick, 112 mm dia

#### Weight

0.7 kg

#### 55X30 Backscatter Section

The Backscatter Section is used both in one-component differential mode and in two-component two-color mode. On the side of the section mounts the 55X34 PM Optics with one or two PM Sections.

#### Backscatter Mirror

First surface type, coated

#### Backscatter Aperture

50 mm

#### Dimensions

70 mm width, 112 mm dia

#### Weight

1 kg

#### 55X32 Beam Translator

The Translator provides continuous variation of beam separation with stops at 13 mm, 26-mm and 39-mm separation for one-channel operation. The same section is used in two channel systems.

#### Beam Translator, Internal Separations

One channel: 13, 26, 39 mm

Two channels: 9.2, 18.4,

27.6 mm

#### Prisms


Reflection prisms

#### Dimensions

35 mm width, 112 mm dia

#### Weight

0.8 kg

#### 55X27 Beam Splitter

The 55X27 Beam Splitter is of a color separating type and splits the incoming two-color beam into a green and a blue parallel to the axis. This module is used in two-color setups.

#### Beam Splitter

Color sensitive prism type, 488/515 nm

#### Beam Separation

Green and blue beams separated 30 mm from optical axis in the same plane

#### Dimensions

35 mm width, 112 mm dia

#### Weight

0.7 kg

#### 55X29 Bragg Cell Section

This section will optically shift the incoming coaxial beam 40 MHz relative to the other. Both in case of one-channel and two-channel measurements, frequency shift can be added to the system by installing one 55X29 Bragg Cell Section. The section has been designed so that the optical path lengths of the two beams are equal.

#### Bragg Cell Frequency

40 MHz

#### Signals for Bragg Cell Driver

40 MHz, 1 Vrms at 50  $\Omega$

DC-voltage, 12 V, 300 mA at one BNC input socket

#### Bragg Cell Efficiency


Max. diffracted light 85%

#### Dimensions

70 mm width, 112 mm dia

#### Weight

1.3 kg

#### 55X31 Pinhole Section

The Pinhole Section is used in backscatter measurements (one- or two-channel) and depicts the measuring volume on the pinhole, thereby constituting a very efficient spatial filter eliminating undesired reflections from your flow setup.

#### Optical System

Two achromatic lenses, aperture 50 mm

#### Pinhole Diameter

0.5 mm

#### Pinhole Adjustment Range


$\pm 2$  mm

#### Dimensions

245 length, 112 mm dia

#### Weight

1.8 kg

#### 55X33 Lens Mounting Ring

The ring is used for mounting one of the 55X51-59 Front Lenses. When incorporating the 55X12 Beam Expander in the system this is mounted between the 55X33 and the front lens.

#### Internal Thread

M105x1

#### Dimensions

20 mm thick, 112 mm dia

#### Weight

0.2 kg

# 5X Modular LDA Optics

## 55X34 PM Optics

The photomultiplier optics is used both in forward scatter and backscatter modes. Its objective focusses on the measuring volume (with the aid of an eyepiece).

### Objective

Focal length: 150 mm, aperture 47 mm, achromatic

### Variable Focus

600 mm  $\rightarrow$   $\infty$

### Close-up Lenses

Focal Lengths: 80, 150, 300 mm

### Reference Beam Adapter

Pinhole Diameter

0.1 mm

### Pinhole Adjustment


$\pm$ 2 mm

### Dimensions

122 mm length, 60 mm dia

### Weight

0.5 kg

## 55X36/37/38 Interference Filter

These blue, green, and red narrow band filters are used for increasing the signal-to-noise ratio of the Doppler signals by filtering away unwanted light.

### Wavelengths

55X36: 468 nm

55X37: 514.5 nm

55X38: 633 nm

### Transmission

>75%

### Bandwidth

10 nm

### Effective Filter Diameter

12 mm

### Dimensions

33 mm dia, thread M30x1

## 55X41/42 Mounting Bench

The bench mounts the laser head and the 55X Optics to a rigid system. The bench features fine adjustment facilities for the optical axis and has three adjustable feet. It includes mounting parts for the laser in question.

### Dimension (HWL)

55X41: 60x120x1750 mm

(for Argon-Ion laser)

55X42: 60x120x1400 mm

(for HeNe laser)


### Weight

55X41: 9.0 kg

55X42: 7.5 kg

## 55X43 Tripod

The tripod is useful for mounting the photomultiplier in forward scatter mode. The tripod comprises a spherical joint and holder for two PM Optics to be used in two component reference beam mode.

## 55L10 Photomultiplier

This photomultiplier with extra ordinary sensitivity covers applications with low and medium velocity. It is equipped with an eyepiece and an adjustable pinhole.

### Type of PM Tube

EMI 9658B

### Voltage Range

500 - 1400 V, negative

### Overall Sensitivity

Typical 1450 A/1m

### Upper Frequency

17 MHz at 1400 V

### Objective

Focal length 105 mm, 1:4

### Close-up lenses Focal Lengths

200 mm, 333 mm, 600 mm

### Dimensions

400 mm length, 75 mm dia

### Weight

2 kg

## 55X35 Color Separator

The separator is used in two-color systems and mounts two 55X08 PM Optics.

### Color Separator

Color separating beamsplitter prism

### Separating Wavelengths

488 or 515 nm, efficiency


>98%

### Dimensions (HWD)

37x47x57 mm

### Weight

0.2 kg

## 55X40 Accessories

This set of accessories comprises two supports, a surface mirror, a neutral density filter, and a test objective. The mirror is valuable when checking the optical alignment of the system, and the test objective is used for checking the beam intersection in measuring volume.

### First Surface Mirror

Diameter 100 mm, adjustable optical axis ca.  $\pm$ 10 mrad 2 axes

### Neutral Density Filter

Attenuation: x 100 (488 to 514.5 nm)

### Test Objective

Focal length 4 mm

### Weight

Whole set 2.0 kg

## 55X44 Instrument Box

The instrument box has room for the various 55X modules and is of a flexible design.

### Finish


Mahogany

### Dimensions (HWD)

210x420x390 mm

### Weight

9 kg

## 55X51-59 Front Lens

A complete line of planoconvex and achromatic front lenses are available. The focal length corresponds to the distance to measuring volume (in air)

### Coating


Hard AR-coating, < 0.3% per surface

### Diameter

115 mm, thread M105x1


Type No.	Focal Length mm	Aperture mm	Weight kg
55X51, Planoconvex	300	94	0.2
55X52, Planoconvex	600	94	0.2
55X53, Planoconvex	1000	94	0.2
55X54, Planoconvex	1500	94	0.15
55X55, achrom.	80	48	0.6
55X56, achrom.	160	48	0.4
55X57, achrom.	310	80	0.5
55X58, achrom.	600	94	0.6
55X59, achrom.	1200	94	0.5

# LDA Signal Processing

**55L11 Diode Detector**  
The 55L11 Diode Detector comprises an aperture, a beam splitter, two photo diodes and a circuit board enclosed in a housing. The complete assembly is mounted in a support which allows one complete 360° rotation of the Diode Detector.

**Sensitivity**  
Typical 2V/mW at 200 kHz input  
max. 400 mVpp  
**Bandwidth**  
2 kHz to 50 MHz (3 dB)  
**Beam Splitter**  
Anti Reflex coated, polarized for 633 nm  
**Power Supply Voltage**  
18 V DC, 40 mA  
**Output Impedance**  
50 Ω  
**Length**      **Diameter**  
125 mm      31.5 mm  
**Weight**  
0.3 kg

**55L17 Seeding Generator**  
The operation of a laser Doppler Anemometer depends on the availability of particles to scatter incident light on to a photo-detector. The particles must have a certain size to produce good Doppler signals and yet be small enough to follow the turbulence of the medium.

The DISA Type 55L17 Seeding Generator constitutes a valuable aid when performing LDA measurements in gaseous flows. The Seeding Generator comprises a carefully developed liquid atomizer and separator. The liquid to be atomized is poured into the Seeding Generator through a filler hole. Compressed air (gas) is fed through a reduction valve (not delivered by DISA) to the IN hose connector. A highly seeded gas flow for your investigation setup is available at the OUT hose connector.

**Particle Density**  
Approx.  $5 \times 10^6$  particles/cm<sup>3</sup>  
**Working Pressure** (= differential Pressure between IN and OUT hose)  
0 to 150 kPascal (0 to 1.5 bars)

**Maximum Supply Pressure**  
600 kPascal (6 bars)  
**Testing Pressure**  
1.5 MPascal (15 bars)  
**Nozzle and Separator Material**  
Glass  
**Bowl Material**  
Transparent polycarbonate plastic  
**Atomizing Liquid Volume**  
Max. 110 cm<sup>3</sup>  
**Hose Connectors**  
8-mm dia.  
**Dimensions (HWD)**  
170x160x90 mm  
**Weight**  
1.7 kg

**55L67 Optical Unit**  
The type 55L67 Optical Unit meets the need for a completely integrated optical unit that is easy to set up and requires no adjustment when beginning another measurement.  
**Beam Separation**  
45 mm or 22.5 mm  
**Lens Focal Length**  
120 mm and 300 mm  
**Optical Surfaces**  
Hard coated anti-reflex front lens and beam splitter  
Quartz coated metal surface mirrors  
**Lead**  
UNF 1" x 32 (laser head standard)  
**Dimensions**  
60 mm dia. x 60 mm  
**Weight**  
1.45 kg  
**Polarization Rotator**  
Code No. 9058C610

**LDA Signal Processing**  
The instantaneous velocity  $V_x$  is represented by the Doppler heterodyning frequency  $f_D$ :

$$f_D = V_x \cdot \frac{\lambda_e}{2 \sin \theta / 2}$$

In all cases, however, demodulation and processing of the Doppler signal must be adapted to the type of flow under investigation. Flows having a low density of scattering particles will typically produce a signal consisting of short bursts with a good S/N ratio, separated by no-signal intervals. The LDA Counter Processor is ideal for processing this signal: It filters


the photodetector signal to remove the DC component (DC pedestal) and the high-frequency noise, and thereafter uses the filtered Doppler bursts to gate a high-frequency clock during a known number of Doppler periods,  $n$ . After the gate has closed, the period time  $T_D$  of the Doppler signal can be calculated by dividing the content of the clock register  $m_c$  by  $n$ . The Doppler frequency  $f_D$  can thus be found from

$$T_D = \frac{m_c \cdot T_c}{n} \Rightarrow$$

$$f_D = \frac{1}{T_D} = \frac{n}{m_c \cdot T_D}$$

$T_c$  being the period time of the clock (Fig. 1).

Fig. 1. Counting Technique


The use of high-speed logics and a high-frequency clock makes it possible to use LDA Counter Processor to measure flow velocities far into the supersonic range.


In flows having a high density of scattering particles it is possible to obtain Doppler signals in which coherent bursts are only separated by brief dropouts with random phase changes. The semi-continuous nature of the signal assures that no sudden changes in frequency will occur between consecutive bursts and allows the use of a tracking-type signal processor, the Doppler Frequency Tracker. The Doppler

Frequency Tracker continuously compares the filtered Doppler signal with the output of a local oscillator, and any frequency difference between the two signals is used to lock the local oscillator to the Doppler frequency. Accordingly, by monitoring the continuous signal from the local oscillator it is possible to make real-time flow velocity measurements (Fig. 2). Since the Doppler Frequency Tracker acts as a narrow-band tracking filter it is capable of processing even noisy photodetector signals that may occur under adverse measuring conditions.

Fig. 2. Frequency Tracker

# 55N20 Doppler Frequency Tracker

The DISA Type 55N20 Doppler Frequency Tracker

## Advantages

- Total tracking frequency from 1 kHz to 10 MHz
- Covered by a total of seven separate overlapping ranges
- Each range has a 1 to 10 frequency ratio
- LED indication of tracker range selected
- Phase-locked-loop controlled VCO
- VCO frequency equal to Doppler frequency
- Lock detector with automatic lock-in search
- Spectrum analyzer function for simplified Doppler frequency detection
- Computer operated remote tracking range and filter control
- Input selector for choice of photomultiplier or diode detector
- Digital output panel for buffer and computer-interface
- Analog output with selectable low-pass filters
- Features a microprocessor for mean velocity calculation
- Direct m/s velocity display
- Calibration factor selection capability
- Averaging time selection of mean velocity display
- Frequency shift compensation facility

## Introduction

The DISA Type 55N20 Frequency Tracker is a completely new design in laser Doppler anemometry (LDA) instruments. It measures the instantaneous fluid flow velocity and presents direct m/s display. An analog output voltage is available for signal processing and digital outputs are available for computer analysis, as well. The 55N20 Frequency Tracker is extremely simple to operate and features a logical front-panel layout. The rear panel contains all connecting terminals to provide versatile operation of the 55N20 Tracker with other instruments such as a frequency shifter or a mini-computer.

## Modular Design

The 55N20 Frequency Tracker is made up of the following modules:

- 55N21 Main Frame
- 55N22 Power Supply
- 55N23 Tracker
- 55N24 Display Module (optional)
- 55L97 High Voltage Supply (optional)

The 55N24 Display Module calculates the mean velocity while the 55L97 HV Supply provides voltage for photomultiplier operation.


## Principle

The incoming Doppler signals are amplified and filtered. These signals, together with the VCO-out signal, are fed to a phase detector. The phase detector signal goes to the Loop Integrator which in turn controls the VCO. The Analog Output signal originates at the Frequency-to-Voltage Converter and passes through the Output Filter. The VCO-control signal is also present here. The VCO output to the Digital Out circuit is formed by the counting circuit which counts the VCO over one complete time-period (both, positive and negative level shifts of the VCO are counted). In the display Module the VCO counts are converted to the Mean Velocity Out display, with Frequency Shift compensation. The Lock Detection circuit detects whether the Phase Locked Loop is in or out of lock. If out of lock, the Search function will startup and search for lock-in, to begin with close to the fre-

quency where the lock was missed, and then with increasing cycle amplitude to cover the entire tracking range. The Lock Detection is delayed 500 VCO signal periods before searching begins. When the search is concluded the phase locked loop goes directly into lock. The 55N20 Tracker can function as a Spectrum Analyzer; i.e. in the SCAN mode the VCO is swept through the range and the Lock Detector circuit is connected to the Mean Square Calculation. When displayed on an oscilloscope the Spectrum Out signal will show the actual Doppler spectrum.

## Technical Data

- 55L97 High Voltage Supply
  - Voltage Range 200 to 2000 V, negative polarity
  - Ripple Max 35 mVrms
  - Voltage Meter Range 0 to 2 kV
  - Anode Current Meter Range 0 to 100  $\mu$ A

Simplified block diagram of the 55N20 Frequency Tracker

### 55N23 Tracker

Photodetector Input  
Impedance: 50  $\Omega$   
Diode Detector: +18 V DC at 40 mA ( $R_f = 150 \Omega$ )

Preamplifier Gain  
40 to 90 dB  
Bivalent Self-Noise  
Typical 1.2 nVA/Hz, referred to input

Electronically Switched Filters  
Controlled by RANGE selector

Monitor Out Terminal  
50  $\Omega$  out  
Use 50  $\Omega$  cable and high impedance load||max. 100 pF

Max. output level 400 mVpp

Type of Tracker Circuit  
Phase-locked-loop

Lock Detection  
Automatic search function

Frequency Tracking Ranges  
1 to 10 kHz, 3 to 33 kHz,  
10 to 100 kHz, 33 to 333 kHz,  
0.1 to 1 MHz, 0.3 to 3.3 MHz,  
and 1 to 10 MHz

Upper Frequency for Flow Fluctuations

Maximum small signal tracking frequencies in resp. ranges:  
0.74 kHz, 2.4 kHz, 7.4 kHz,  
24 kHz, 74 kHz, 240 kHz,  
0.74 MHz

Capture Range in resp. Freq. Ranges  
 $\pm 0.74$  kHz,  $\pm 2.4$  kHz,  $\pm 7.4$  kHz,  
 $\pm 24$  kHz,  $\pm 74$  kHz,  $\pm 240$  kHz,  
 $\pm 0.74$  kHz

Maximum Slew-Rate in Resp. Ranges  
210 kHz/sec, 630 kHz/sec,  
6.3 kHz/msec, 63 kHz/msec,  
630 kHz/msec, 6.3 MHz/msec,  
63 MHz/msec

VCO Frequency Ranges  
Equal to tracking ranges

Analog Output  
Voltage range: 1 to 10 V for all seven ranges  
Accuracy: 1% of full range  
Impedance: 100  $\Omega$

Output Filter  
First order low-pass filter.  
Four settings in each tracking range. In 10 MHz range:  
1.9 kHz, 5.9 kHz, 19.4 kHz and  $\infty$

Spectrum Output  
Voltage range: 0 to 10 V  
Impedance: low

VCO Control Output  
Voltage range: 1 to 10 V for all seven ranges  
Accuracy: 4% of full range  
Impedance: low

Sweep Input  
Voltage range: 1 to 10 V  
Impedance: high

Lock Detector signal  
Open collector with 2 k $\Omega$  pull-up  
Digital Input Terminal  
Remote setting of range and output filter  
Frequency shift information  
Digital Output Terminal  
Setting of range and output filter (values may be a result of local or remote settings)  
Doppler frequency - 8 bits  
Frequency shift information  
Computer handshake signals

### 55N24 Display Module

Mean Velocity Display  
3 digits with automatic selection of sign, decimal point after first digit  
Power of 10 display: range from -9 to +9  
Calibration Factor Range  
 $10^{-3}$  to  $10^2$  ms $^{-1}$ /MHz  
Integration Time Range  
1 to  $9 \times 10^2$  sec  
Range of Frequency Shift Input  
0 to  $\pm 9$  MHz  
Mean Velocity Output Terminal  
Average Doppler frequency - 16 bits  
Computer handshake signals


### 55N22 Power Supply

Line Voltage  
For operation on any nominal line voltage from 100 to 129 and 200 to 259 V AC, 50 to 60 Hz.  
Voltage selection by means of selector switch  
Power Consumption  
60 VA  
Fuse  
0.5/1.25 A at 220/110 V

### Additional Data

Dimensions (HWD)  
133 x 440 x 420 mm  
Weight 11.6 kg

# 5L90a LDA Counter Processor


The DISA Type 55L90a LDA Counter Processor

## Features

- Digital display of mean velocity in m/sec
- Auto ranging over full Doppler frequency from 2 kHz to 100 MHz
- Four modes of operation: Fixed, combined, variable and transit
- Advanced group of data validation circuits: amplitude domain validation with three Schmitt triggers
- Time domain validation by odd number (5/8) fringe count
- Automatic selection of number of fringes between 2 and 256
- True 500 MHz clock
- Total burst mode
- Digital display of data rate and percent validated data
- Large particle rejection circuit
- Data interface and buffer provided by 57G20 Buffer Interface
- High-pass and low-pass filters
- Digital and analog outputs

## Introduction

Measuring the velocity of a fluid flow by means of a counting technique is the latest development in laser anemometry. Fundamentally, the Counter times a scattering particle over a known number of interference fringes within the measuring volume, thus enabling evaluation of instantaneous velocity.

Doppler signal from a signal scattering particle and the output from the Schmitt triggers.

## Principle

The separation between fringes is known from the geometry of the setup and wavelength of the laser light. A typical Doppler signal for a scattering particle crossing the center of the measuring volume is shown in Fig. 1. This signal is band-pass filtered to remove high-frequency and low-frequency noise, amplified and allowed to trigger a Schmitt trigger. Fig. 3 shows the resulting output from the Schmitt trigger. The period of this square wave is equal to the time it takes for a scattering particle to cross one fringe.

If the LDA Counter is operated in mode 1 (fixed  $N_f$  mode used mainly for almost continuous Doppler signals), operation is as follows: The Schmitt trigger output is fed to a control circuit which enables two high-speed counters (High Count and Low Count both counting a 500 MHz clock frequency) on the same zero crossing. After five fringe counts have been completed at the Fringe Counter input the Low Count is disabled. After eight fringe counts have occurred at the Fringe Counter input the High Count is disabled. A Comparator circuit then compares the Low Count with 5/8 High Count. An "accept" or "reject" command is released according to the desired accuracy (set by COMP. ACCURACY selector).

In certain cases false data can pass through the Comparator. A Two Level Validation circuit comprising two Schmitt trigger levels and a Sequence Detector eliminates such false validation. If the comparison results in an "accept" command, the High Count, being directly proportional to the time it takes for the scattering particle to cross eight fringes, is inverted giving an output which is directly proportional to the instantaneous velocity. This result is passed to further processing circuits such as the 55L94 Mean Velocity Computer and 55L93 D/A Converter.

The LDA Counter is operated in mode 3 (variable  $N_f$  mode used mainly for "single burst" Doppler signals), operation is as follows:


The High Count register also in this mode is enabled on the zero crossing of the Schmitt trigger. It stays enabled as long as the Doppler signal is above a certain level. In this way the "burst time" is measured and at the same time the number of periods is counted by the Fringe Counter. These two data allow an interfaced computer to calculate the velocity and also to make very accurate corrections for bias. The advantages of the two above-mentioned operating modes are combined in mode 2: Combined mode.

It should be mentioned that in a mode of operation, the transit mode (mode 4), is used for measuring the velocity of other particle transit phenomena. The time from one trigger event to the next (two channels) is measured and the inverted time is output.

#### Technical Data

The 55L90a LDA Counter Processor consists of the following modules which may be ordered separately, if desired:

- 55L96 Counter Module
- 55L97 High Voltage Supply Module
- 55L91 Data Rate Module
- 55L93 D/A Converter Module
- 55L94 Mean Velocity Computer Module
- 55L99 Power Supply

Block diagram showing principle of operation.

#### 55L96 LDA Counter Module

##### Modes of operation

1. Fixed  $N_f$  mode
2. Combined mode
3. Variable  $N_f$  mode
4. Transit mode

Input Impedance  
50  $\Omega$

Input Frequency Range  
2 kHz to 100 MHz

Calibration Factor Range  
1.3 to 44 m/sec/MHz

(55X LDA Optics)

Velocity Range  
3 mm/sec to 2000 m/sec

(approximate values for DISA LDA optics - dependent on geometry of setup)

Burst Time/Transit Time Range  
(mode 3, 4)

64 nsec to 4.2 msec

Number of Fringes  
2 to 256

Low Pass Filter Settings  
256 kHz

1, 4, 16, 100 MHz (3 dB)

Roll-off  
60 dB/decade

High Pass Filter Settings  
1, 4, 16, 64, 256 kHz

1, 4, 16 MHz (3 dB)

Roll-off  
40 dB/decade

Threshold Window  
20 dB in 1 dB steps

Instrument Limited Sample Rate  
1.0  $\mu$ sec + measuring time

#### Accuracy ( $\approx$ SNR)(fixed $N_f$ )

1 % at 40 MHz

2.5 % at 100 MHz

#### Digital High Pass Filter

Selects a minimum frequency or maximum measuring time

#### 55L97 High Voltage Supply Module

Voltage Range  
200 to 2000 V negative polarity

Ripple  
Max. 35 mV rms

Voltage Meter Range  
0 to 2 kV

Anode Current Meter Range  
0 to 100  $\mu$ A

#### 55L94 Mean Velocity Computer Module

##### Digital Display

3 digits of velocity  
1 digit of power of 10 with automatic selection of sign

Range of Scale Factor Setting  
0 to 9999, Scale setting for velocity display in MHz: 0150

#### 55L91 Data Rate Module

Data Rate Display  
.000 to 999 kHz

Validated Display  
000 to 999 parts per thousand (based on last 1000 samples)

#### 55L93 D/A Converter Module

Output Resolution  
0.2 % of selected range

Output Impedance  
600  $\Omega$  (DC)

Low Pass Filter  
300 kHz (3 dB)

Roll-off  
20 dB/decade

Output Range  
0 to 1.25 V

Output Calibration Factor

$$f_0 \text{ (MHz)} = \frac{100}{2A} \times V_{out} \text{ (V)}$$

A = front panel gain setting

#### 55L99 Power Supply

Line Voltage  
For operation on any nominal line voltage from 100 to 129 V and 200 to 259 V AC, 50 to 60 Hz

Voltage selection by means of selector switch


Power Consumption  
180 VA

#### Additional Data

Dimensions (HWD)  
190 x 440 x 420 mm

Weight  
16 kg

# 55N10 LDA Frequency Shifter

The DISA Type 55N10 LDA Frequency Shifter.

## Features

- \* Modular design — 1- or 2-channel version
- \* Directly suitable for operation with DISA counters/trackers
- \* Frequency shift range  $\pm 9$  MHz
- \* Provides 40-MHz signal for Bragg cell
- \* Binary outputs facilitate computer interfacing

## Applications

Optical and electronic frequency shifting allows measurements of highly turbulent and reversing flows such as:

- \* Flow profiles with velocity reversals
- \* Oscillatory flows of zero mean velocity
- \* Vortex formation in wake flow
- \* Natural convection flow
- \* Unstable boundary layer flow

The DISA 55N10 LDA Frequency Shifter is operated in conjunction with LDA signal processors such as DISA counters/trackers. Because of its modular design both a 1-channel and a 2-channel version are available. Effective frequency shifts from 10 kHz to 9 MHz are provided on both channels, with individual frequency selectors.

The direction of the effective shift can be reversed relative to the sense of flow without re-adjustment or turning of the optics. A 40 MHz signal is provided for driving the Bragg cell for optical frequency shift. Terminals on the rear panel carry binary signals relating to the knob settings for the channels in question, for computer interfacing.

## Technical Data

The 55N10 LDA Frequency Shifter comprises the following modules:

- 55N11 Frame with Mother Board
- 55N12 Mixer Unit
- 55N13 LM Exciter Unit
- 55N14 Frequency Shift Generator Unit
- 55N15 Fixed Divider Board
- 55N16 Power Supply

A one-channel setup comprises the following units: 55N11, 55N12, 55N13, 55N14, 55N15, 55N16. Expansion to a two-channel setup is effected by adding two extra sections: 55N12 and 55N14.

### 55N11 Frame with Mother Board

Consists of two D subminiature 9-spole female sockets with locking springs. The sockets are designated FREQUENCY SHIFT, DIGITAL OUT, CHANNEL 1 and CHANNEL 2. The binary output of the DIGITAL OUT sockets indicates the setting of the FREQUENCY SHIFT GENERATOR UNIT. Output level is standard TTL with a fan-out of 4 per pin.

### 55N12 Mixer Unit

PM Input  
Input impedance: 50  $\Omega$   
Input voltage range: 10  $\mu$  Vrms to 2.5 mVrms  
Input frequency range: 10 MHz to 50 MHz  
Equivalent input white noise spectrum: Typically 4 nV/ $\sqrt{\text{Hz}}$

$f_{LO}$  Input from 55N14  
Input impedance: 500  $\Omega$  / 560 pF  
Input voltage range: 70 to 100 mVrms  
Input frequency range: 30 MHz to 50 MHz

Mixer Out  
Output impedance: 50  $\Omega$   
Conversion gain: Typically 15 dB  
Output frequency range: 1 kHz to 20 MHz  
Output voltage range: 56  $\mu$ V to 14 mV  
Anode current indication signal: Max. 200  $\mu$ A.

### 55N13 LM Exciter Unit

General  
Output frequency  $f_0$ :  
40 MHz  $\pm$  4 kHz  
Frequency drift: 12.5 ppm  
Frequency drift or frequency instability is a term expressing the maximum change of frequency relative to 40 MHz taking place over a temperature range from 0° to 60°C.

40 MHz Output to Bragg Cell  
Output voltage (adjustable):  
0.1 Vrms to 1.7 Vrms into 50  $\Omega$ .  
Output voltage, DC, unloaded:  
15 VDC  
Output voltage at max. 300 mA:  
12 VDC  
(load 50  $\Omega$  plus built-in 10  $\Omega$  limiting resistor)

40 MHz Output to 55N15  
Fixed Divider Board:  
Output Voltage: 40 mVrms into 50  $\Omega$

### 55N14 Frequency Shift Generator

Output frequency: 40 MHz  $\pm$  (1, 2, ..., 9) x (10, 100, 1000) kHz.

Long term frequency drift: As for 55N13 (frequency locked)  
Short-term frequency ripple (jitter): Typically  $\pm$  50 Hz.  
Output voltage: 70 to 100 mVrms (unloaded)  
Output impedance: 50  $\Omega$   
Spurious frequencies: Better than -30 dB relative to  $f_{LO}$ .  
FREQUENCY IN KHz selector: 9 position switch selects factors between 1 and 9 kHz  
MULTIPLIER selector: 6 position switch selects both frequency multiplication factor and sign of frequency shift (+ when  $f_{LO} > 40$  MHz, - when  $f_{LO} < 40$  MHz)

Input frequency: 2.5 kHz (40 MHz divided by 16000)  
Input voltage: TTL level from 55N15

### 55N15 Fixed Divider Board

Input frequency: 40 MHz from 55N13  
Input voltage: 40 Vrms  
Input impedance: 50  $\Omega$   
Output voltage: TTL level to 55N14  
Output frequency: 2.5 kHz (40 MHz divided by 16000)

### 55N16 Power Supply

Line voltage: 100, 120, 140, 200, 220, 240 VAC, 50/60 Hz. Stabilized against line voltage variations of  $\pm$  10%.  
Power consumption: 40 VA  
Output voltages: Stabilized +5 V and +15 V  
Fuse: Slow-blow, 5 x 20 mm; 110 VAC: 0.63 A  
220 VAC: 0.32 A

## Additional Data

Ambient temperature range: +5° to +40°  
Dimensions (HWD): 150 x 440 x 420 mm  
Weight: 8.5 kg (2-channel version).

# 57G20 Buffer Interface

## Features:

- \* Provides communication between DISA Scientific Research Equipment and a PDP-11 Computer
- \* Real-time data logging rate up to 1 MHz (10 MHz optional)
- \* Simultaneous events logging in multichannel setup
- \* Buffer of any length (expandable buffer)
- \* Four double line input channels
- \* Eight outputs for control of computer-operated auxiliary equipment
- \* Ideal for operation with 55L90 Counter and 55N20 Tracker
- \* Highly versatile

## Introduction

The DISA 57G20 Buffer Interface system comprises a frame and a line of circuit boards for interface signal processing.


## Units:

- 57B02 Main Frame
- 57G106 Control Logic
- 57G113 Output Demultiplexer
- 57G120 Input Multiplexer
- 57G127 Output Logic Buffer,  $\frac{1}{2}k$
- 57G134 Buffer Memory,  $\frac{1}{2}k$
- +5V/-12V DC-to-DC Converter
- Coincidence Filter
- Serial Transmission Line

The 57B02 Main Frame has a sink front panel and a rear panel carrying a power supply terminal. The Frame has a 40-channel series-operating backplane and will accommodate up to 10 circuit boards.

The 57G106 Control Logic circuit board controls data reading: It gives a priority order to the input channels in case of simultaneity, enables single or double input on the channels, matches input data to the buffer registers, and controls timing interaction and information in case of buffer overload. A socket for system self-testing is provided.

The 57G113 Output Demultiplexer is a self-contained communications circuit board. It demultiplexes output from the PDP-11 into eight channels, two of which are used to control functions (time base, ranges, etc.) of the DISA Counter/Tracker. The six are available for your computer controlled measuring setup, e.g. traversing.

The 57G20 Buffer Interface comprises a cabinet, flat cables, and circuit boards. Holders for rack mounting are also included.

The 57G120 Input Multiplexer board has four double channels each of which connects through a 40-conductor flat cable to your DISA Counters/Trackers (ordering no.: 06A214 Flat Cable).

The 57G127 Output Logic Buffer board comprises a variable length buffer register (maximum  $\frac{1}{2}k$  16-bit words), control logic for output data, and a +5V/-12V DC-to-DC converter. This board buffers random asynchronous signals from your measuring setup into synchronous signals matching the Computer.

The 57G134 Buffer Memory boards are used to expand your buffer capacity, each board by  $\frac{1}{2}k$ .

The +5V/-12V DC-to-DC Converter board powers the 57G134 Buffer Memories.

The Coincidence Filter board is intended mainly for two- or three-dimensional velocity vector measurements. It blocks the Input Multiplexer (except in case of data coincidence from the actual channels). This feature is of importance in cross-correlation computations. The coincidence time window is computer controlled.


The Serial Transmission Line (medium connection) board enables operation of your measuring setup in conjunction with a remote computer, e.g. a central computer.

## Technical Data

- 57B02 Main Frame Backplane
  - 2x41-pole board connectors
  - Number of Circuit Boards Max. 10
  - Power Supply +5V, max. 10A
  - Rear panel socket: Fischer F-400A002+E
  - Dimensions (HWD) 133 x 440 x 420 mm
  - Weight 4 kg
- 57G106 Control Logic Control Functions
  - Computer handshake, priority order to input channels, single or double input to channels, adaptation of input data to buffer register, overflow indication of buffer, system self-testing.
  - Power Supply +5V, 250mA
- 57G113 Output Demultiplexer Input Terminal
  - 40-pole 3M 3432-1002 (mating plug type: 3417-3000 - 3490-5)
- Demultiplexer Outputs
  - 6 TTL compatible external outputs (16 pin DIL sockets)
  - 2 outputs employed inside instrument
- Output Signal
  - 12 bit word + reference ground and new data ready line
- Power Supply +5V, 400mA

- 57G127 Output Logic Buffer Buffer Register
  - 16-bit words, maximum length  $\frac{1}{2}k$
  - Input Rate Max. 1 MHz (10 MHz optional)
  - Control Logic Generates handshake to PDP-11
  - DC-to-DC Converter Output +5V to -12V, 100mA (1 MHz version only)
  - Output Terminal 40-pole 3M 3432-1002 (mating plug type: 3417-3000 + 3490-5)
  - Power Consumption +5V, 1.3 A
- 57G134 Buffer Memory Buffer Register
  - 16-bit words, maximum depth  $\frac{1}{2}k$
  - Input Data Rate Max. 1 MHz (10 MHz optional)
  - Power Consumption +5V, 500 mA
  - 12V, 90mA (1 MHz version only)
  - DC-to-DC Converter Power Supply +5V, depending on load (max. 3.5 A)
  - Output Voltage -12V, max. 1.2 A

## Educational LDA System

### Features

- \* Extremely easy to set up
- \* Very easy to operate
- \* Measured velocity vector rotatable through 90°

This tracker-based LDA System operates in the reference beam mode featuring extreme ease of optical alignment: Mount the 55L67 Optics on the laser head, select the proper front lens, and direct the 55L11 Diode Detector onto the measuring point with the reference beam impinging the Diode aperture. The system is well suited for water flow applications and the optical Doppler pedestal removal secures a very high S/N ratio.

### System Data

**Modes of Operation**  
Reference beam mode

**Flow Medium**  
Water


**Velocity Range**  
2 mm/s to approx. 50 m/s

**Turbulence Intensity**  
Max. 30%

**Distance to Point of Measurement**  
120 mm resp 300 mm

### System Units

One 55N20 Doppler Frequency Tracker  
One 55L67 Optical Unit  
One Retarder Disc 9055C6101  
One 55L11 Diode Detector

## Low Velocity LDA System

### Features

- \* Easy-to-operate system
- \* High sensitivity of photomultiplier
- \* Measured velocity vector rotatable through 90°

The low velocity tracker based system will operate in differential Doppler mode in water flows or moderately seeded air flows. Due to the extremely high sensitivity of the PM tube this system is suited for flows with a particle concentration corresponding to only a single particle in the measuring volume at a time.

### System Data

**Modes of Operation**  
Differential Doppler mode

**Flow Medium**  
Air (gas)  
(Water)

**Velocity Range**  
Approx. 2 mm/s to 50 m/s


**Turbulence Intensity**  
Max. 30%

**Distance to Point of Measurement**  
120 mm resp. 300 mm

### System Units

One 55N20 Doppler Frequency Tracker, incl.  
55L97 High Voltage Supply Unit  
One 55L67 Optical Unit  
One 55L10 Photomultiplier

## Coaxial Laser System

### Features

- \* Compact LDA optics
  - \* Tracker-based or counter-based system
  - \* Rotatable optics through 360°
  - \* Wide range of front lenses
- The coaxial laser system has adapters for mounting of most commercial laser heads. The LDA optics is light and compact, but the limited power of today's coaxial laser restricts the velocity range to low and medium. The flexibility of this system ensures measurements of one and two component forward scatter setups - three different beam separations - nine different front lenses are available. You have a choice of selecting a tracker-based system which is well suited for liquid measurements or a counter-based system optimized for gas flows. In both cases operation of the systems is easy, featuring a mean velocity display direct in m/s.

### System Data

**Modes of Operation**  
One component forward scatter in differential or reference beam mode  
Two component forward scatter reference beam mode

**Flow Medium**  
Tracker-based System: Liquid flow and highly seeded air flow  
Counter-based System: Low seeded gas flow

**Velocity Range** (5 mW coaxial laser)  
Tracker-based System: Approx. 2 mm/s to 50 m/s  
Counter-based System: Approx. 2 mm/s to 100 m/s

**Turbulence Intensity**  
Tracker-based System: Max. 30%  
Counter-based System: Max. 50%

**Distance to point of Measurements**  
80 mm to 1.5 m

### System Units


One or two 55N20 Doppler Frequency Tracker, incl.  
55L97 High Voltage Supply or  
One or two 55L90a LDA Counter Processor  
One 55X51-59 Front Lens  
One 55X23 Support  
One 55X24 Beam Splitter  
One 55X26 Beam Splitter  
One 55X28 Beam Displacer

55X24 + 55X28 for single ch. diff./ref. mode  
55X24 + 55X26 for two ch. ref. mode  
55X26 for single ch. diff. mode

One 55X33 Lens Mounting Ring  
One or two 55X08 PM Section  
One or two 55X34 PM Optics  
One or two 55X38 Interference filter (at high ambient light level)  
One 55X19 Laser Adapter  
One 55X43 Tripod


55N20 Doppler Frequency Tracker

55L90a LDA Counter Processor

# Backscatter LDA System


Equipped one velocity component anemometer your flow setup from only velocity vector through 180° velocity range and intensity range flows with zero mean velocity. The LDA system covers flow studies. The LDA equipped with optical shift section and backscatter mounted together with 15 mW laser constituting a highly stable setup. Investigation setup is from one side only, and will measure the fluid flow with extreme accuracy. Optical and electronic shifting the velocity system has great flexibility. When operating the 55N20 Tracker of the mean flow direction, (indicating sign). Any frequency shift is automatically compensated as both a Doppler system or a counter-based system the frequency shift introduces a few controls only.


Operation  
 Component forward  
 and backscatter in  
 all Doppler mode  
 in  
 liquid  
 range  
 300 m/s to + 300  
 intensity  
 Point of Measure-  
 1.2 m

## System Units

- One 55L90a LDA Counter Processor
- or
- One 55N20 Doppler Frequency Tracker
- One 55N10 Frequency Shifter (1-ch. version)
- One 55X55-59 Achromatic Front Lens
- One 55X08 PM Section
- One 55X12 Beam Expander
- One 55X21 Cover and Retarder
- Two 55X23 Supports
- One 55X24 Beam Splitter, pol. 633
- One 55X28 Beam Displacer
- One 55X29 Bragg Cell Section
- One 55X30 Backscatter Section
- One 55X32 Beam Translator
- One 55X33 Lens Mounting Ring
- One 55X34 PM Optics
- One 55X38 Interference Filter
- One 55X42 Mounting Bench (for 15 mW laser)
- One 55X43 Tripod


*Tracker-based system with Frequency Shifter*

*Counter-based system with Frequency Shifter*

# Two-Color LDA System

## Features

- ★ Two-dimensional flow field investigations
- ★ Unseeded or sparsely seeded fluid flow measurements
- ★ Reynolds shear stress determination (correlation)
- ★ Supersonic flow measurements

Providing simultaneous measurements of two perpendicular velocity components, the Two-Color LDA System operates in differential Doppler backscatter mode. It utilizes the two strongest laser lines of an argon ion laser (green and blue lines) and the high output of the laser is of importance in supersonic measurements. The electronics comprises two 55L90a Counters and a two-channel version of the 55N10 Frequency Shifter.

## System Data

### Modes of Operation

Two component two-color forward scatter and backscatter modes

### Frequency Shift

Optical and Electronic

### Flow Medium

Gases and liquids

### Velocity Range

Approx. -300 to +1500 m/s (Mach 5)

### Turbulence Intensity


Infinity

### Distance to Point of Measurement

80 mm to 1.2 m

## System Units

- Two 55L90a LDA Counter Processors
- One 55N10 Frequency Shifter (2-ch. version)
- Two 55X08 PM Sections
- One 55X12 Beam Expander
- One 55X20 Cover and Retarder
- One 55X22 Beam Waist Adjuster
- Two 55X23 Supports
- One 55X25 Beam Splitter
- One 55X27 Beam Splitter
- One 55X29 Bragg Cell Section
- One 55X30 Backscatter Section
- One 55X31 Pinhole Section
- One 33X32 Beam Translator
- One 55X33 Lens Mounting Ring
- One 55X34 PM Optics
- One 55X35 Color Separator
- One 55X36 Interference Filter (blue)
- One 55X37 Interference Filter (green)
- One 55X41 Mounting Bench (For 5 W Argon Ion Laser)
- One 55X40 Accessories
- One 55X43 Tripod
- Two 55X44 Instrument Boxes
- One 55X55-59 Achromatic Front Lens


*Two 55L90a LDA Counter Processors and a two-channel version of 55N10 Frequency Shifter*

# Computer System

rocessing of measur-  
flow studies can be  
d  
elocity measure-  
e computer can  
ne setup, e.g. travers-

checking programs  
rrect operation

ftware  
omputer System is  
ar and constitutes a  
e instrumentation.  
otics matching your  
ring setup is selected  
ed to a tracker-based  
ased processing  
The DISA 55N20  
quency Tracker and  
LDA Counter Pro-  
igital output ter-  
irect connection to  
uffer Interface. The  
face unit can be  
th the desired buffer  
ther control func-  
om this unit the bus  
nected to a PDP  
uter. The configura-  
omputer for your  
requirement can  
ed without any  
nules. For  
measuring data and  
ppy disc module  
orated. For check-  
grams, stored data,  
d data you can  
stance a complete  
ard terminals,  
monitors.

Example of a complete LDA instrumentation comprising Two-Color LDA Optics, two LDA Counter Processors, Frequency Shifter, Buffer Interface and Computer with keyboard terminal.

## Software system

The DISA software system is designed as a number of libraries and as stand-alone programs. The distribution media are punched paper-tape and floppy disc. Having access to an operating system, the user can generate his own calculation programs and controlling facilities.

The following libraries are available:

### Group 0:

This is a collection of test programs, which makes it possible to make hardware diagnostics in a very simple way. This library is always provided as a collection of stand-alone programs.

TEST: starts collection of data, and checks different facilities.

HELPME: is a separation test program. Using this program the hardware error area can be determined.

OUT: checks all optional outputs from the interface. Further there is a command system for manual testing.

### Group I:

This is the data processing subroutine library. The primary design criterion has been data-collection speed. The second criterion is the best possible use of the memory area. COLLEC: takes the arguments from the main section or any appropriate subroutine to set up the mode of operation, number of samples, and so on. In the standard version COLLEC supports 2 double channels. COLLEC is device independent. SORT: is an optional subroutine used to sort random data, if more than one channel is used. This procedure is time-consuming, and should be avoided in the case of non-random data.

TRANSF: is a level translating routine used to establish coincidence between the internal compressed format and the user format.

OUTP1: is a data transmission program used to organize the transmission of collected data to a storage device or a background computer.

SCAN: gives high-level control over optional outputs, traversing control, etc.

### Group II:

This is a collection of calculation subroutines. They are all generated in FORTRAN, and they can be supplied as user written routines or changes. They are all called from a main section, and they must have access to group I or a similar library.

DATA: is the high level data transmission program. It calls TRANSF, and calculates data in user units.

DATAI: is the data-list program. It is used as a diagnostic program to check the setup.


MOMEN: calculates the moments of the distribution function.

HISTO: generates the histograms for all available channels.


CORR: calculates possible components of the correlation functions, both auto- and cross-correlations are available.

SPEC: calculates the spectra based on the correlation functions.

computer operated  
ng systems


Print-out on monitor screen of  
the processing program

Data list showing velocity and  
time information

Histogram showing velocity  
distribution

We reserve the

We reserve the right to make, without notice, such changes in our published data as we may deem necessary or desirable.

# DISA

DISA ELEKTRONIK A/S . DK-2740 SKOVLUNDE . DENMARK

Telephone: + 45 2 84 22 11

Cables: DISAELEKTRONIK, Copenhagen

Telex: 35349 disae dk